

PKA™ Articulating Ankle Foot Orthosis

252SKG-

PKA™ APU® Adult Kodel Gun Metal

*Positioning and Walking Ankle Foot Brace
for Custom KAFO or AFO Adaption*

Features-

- › Adjustability in all planes of motion, depending on design chosen
- › Clinically proven for effective ambulation*
- › Compatible with all liner variations that we offer if applicable (sold seperately)
- › Proven superior integrity of the superstructure**
- › Eliminate posterior heel contact
- › Control dorsi-plantar flexion
- › Controls external or internal rotation of the hip when supine
- › Accommodates those at risk of skin breakdown

*** Heel strap positioned inside
the heel connecting bar for
display purposes only*

-Available in Adult, Pediatric, and Infant

The PKA™ Orthosis should only be used in conjunction with custom made or pre-fabricated Knee Orthoses or Knee Ankle Foot Orthoses. The PKA™ Orthosis should not be used as a standalone Ankle Foot Orthosis.

Uses-

- › For those patients who wear their KAFO 24 hours a day, or spend the majority of their time at rest or in bed.



252PKG
-Adult Pad and Strap



250SKT
-Adult Kodel Teal

PKA™ Orthosis is protected under U.S Patent No. 5,593,383

coding help PDAC Approved L-4396 and/or L-1971

**ANATOMICAL
CONCEPTS, INC.**

